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Clinical Benefits: Evaluation Accuracy comparable to experts: MRMC Study

evaluation, and biopsy would not be recommended for an additional 10.8% of benign nodules. The results demonstrated that applying a CAD system would improve clinicians' interpretations and lessen the variability in diagnosis.




Figure 3. Average AUC Without and With CAD of junior group and senior group. Blue = junior group without CAD; Orange = junior group with CAD; Gray = senior group without CAD; Yellow = senior group with CAD.

Reference: Multi-Reader Multi-Case Study for Performance Evaluation of High-Risk Thyroid Ultrasound with Computer-Aided Detection *Cancers* 2020, 12, 373; doi:10.3390/cancers12020373

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Clinical Benefits: Accuracy comparable to experts: Diagnostic Performance

Reference: Reverter JL, Vázquez F, Puig-Domingo M. (2019). Diagnostic Performance Evaluation of a Computer-Assisted Imaging Analysis System for Ultrasound Risk Stratification of Thyroid Nodules. American Journal of Roentgenology, 213: 169-174.

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510(K) Summary: K192903

In accordance with 21 CFR 807.92 the following summary of information is provided:

1. Date Prepared – October 11, 2019
2. Manufacturer
SAMSUNG MEDISON CO., LTD.
3366, Hanseo-ro, Nam-myeon,
Hongcheon-gun, Gangwon-do 25108,
REPUBLIC OF KOREA
3. Primary Contact Person
Ji Yea Lee
Regulatory Affairs Specialist
Phone: +82.2.2194.1594
Fax: +82.31.8017.9576
Email: jiyea722.lee@samsungmedison.com

Secondary Contact Person
Ninad Gujar
Director of Regulatory & Quality
Phone: +1.978.564.8632
Fax: +1.978.564.8677
Email: ngujar@neurologica.com
4. Proposed Device
Common Name: Diagnostic Ultrasound System and Accessories
Trade/Device Name: RS85 Diagnostic Ultrasound System
Additional Marketing Name: RS85 Prestige Diagnostic Ultrasound System
Regulation Name: Ultrasound pulsed Doppler imaging system
Panel/ Regulatory Class: Radiology / II
Product Code: IYN, IYO, ITX
Regulation No./Product Code: 892.1550, 892.1560, 892.1570
5. Predicate Device
 - Predicates
 - RS85 Diagnostic Ultrasound System (K191115) ; Primary
 - HERA W10 Diagnostic Ultrasound System (K182595)
 - HS50/HS60 Diagnostic Ultrasound System (K181336)
 - Reference
 - AIXPLORER & AIXPLORER Ultimate Ultrasound Diagnostic Systems (K173021)
 - Fibroscan (K150949)
 - AmCAD-US (K162574)

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
Medical Expo: 4/15-7/15

[Click here to see your proposal](#)
(you will be able to review it before confirming)

MedicalExpo -- Membership Offer Starting Package

[Fewer details](#)

Billing upon acceptance

Account Manager: Antoine VERHAEGHE

Telephone: +33.4.86.57.42.38

Start date: Sep 1, 2020

Termination date: Sep 1, 2021

Validity termination date: Jul 15, 2020

ME - AmCad BioMed

Item	Geographic Area:	Quantity	Price
ONLINE STAND - 12 MONTHS	Worldwide	1	3,996.00 EUR

HIGHLIGHT - Newsletter

complimentary

Payment due date: 30 days of receipt of invoice

1Y+1.5m Service Fee:

Total number of items:	3,996.00 EUR
10% Discount	-399.60 EUR
Price before tax:	3,596.40 EUR

MedicalExpo Transitions from an Online Exhibition to a B2B Marketplace

By Celia Sampil, Camille Rustici and Erin Tallman

As companies continue expanding the online market by creating their own marketplaces, the year 2020 might be the cut-off for those interested in entering the competition. That's what Vincent Gerard, CEO of the French company VirtualExpo Group which powers MedicalExpo, thinks. Today, the company is officially launching its new online marketplaces.

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
BD&L : Business Model

Licensing

License technology to equipment makers and receive licensing fees or royalty.

Distribution

Leasing

Revenue
Sharing

Outright
Sales

Subscription

Exclusive or regional rights, single & multi-product deals.

Collaborative Marketing

Collaborate with PACS/AI platforms to market and expand access.

2020 INTERNATIONAL FORUM ON NEW TECHNOLOGIES AND CLINICAL DIAGNOSIS OF THYROID ULTRASOUND

13TH NOV, 2020 (1.5 HOURS) LIVE WEBINAR

10:00 - 11:30 A.M. Spain

16:00 - 17:30 P.M. Indonesia | Vietnam

17:00 - 18:30 P.M. Taiwan | Singapore | Malaysia

Taiwan Society of Ultrasound in Medicine

AmCad BioMed

Register to join webinar

Date: Friday, November 13, 2020

Time: 17:00 -18:30 p.m. (CST/ Taipei time)

Speakers

Prof. Yi-Hong Chou

- Vice President of the World Federation for Ultrasound in Medicine and Biology (WFUMB)
- Chair and Professor of Radiology at Yuanpei University of Medical Technology and National Yang Ming University School of Medicine

Agenda

- Opening Remarks
- Special Messages from Prof. Yi-Hong Chou
- The Trend and New technologies for Thyroid Cancer Diagnosis
By Prof. Tsung-Lin Yang MD, PhD
- Artificial Intelligence (AI) and Its Applications to Sonographic Diagnostics of Thyroid Nodules
By Prof. Argon Chen
- The Role of Computer-Aided Diagnosis (CAD) in the Diagnosis and Characterization of Thyroid Nodules
By Prof. Jordi L. Reverter
- Q&A Panel Discussion
- Closing Remarks

Prof. Tsung-Lin Yang MD, PhD

- Professor and Chief, Head and Neck Surgery, Department of Otolaryngology & Head Neck Surgery, National Taiwan University Hospital
- Chief Executive Officer, Research Center for Development Biology and Regenerative Medicine, National Taiwan University

Prof. Argon Chen

- Professor, Industrial Engineering and Mechanical Engineering, National Taiwan University
- Chief R&D Officer, AmCad BioMed

Prof. Jordi L. Reverter

- Senior specialist Endocrinology, Germans Trias i Pujol University Hospital
- Professor of Medicine, Autonomous University of Barcelona, Spain

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