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


Patent: Ultrasound Scatterer Structure Visualization & Analysis


Patent: Tumor Contour Retrieval

Patent: Multi-Layer Classifier & Echogenicity/ Tumor Feature/ Echo Texture Quantification

Patent: Method for Head & Neck Assessment or Intervention


Other Patent Titles	Countries
用於頭頸部評估或介入治療的定位裝置/ POSITIONING APPARATUS FOR HEAD AND NECK ASSESSMENT OR INTERVENTION	China EU (Germany) Japan Taiwan USA
呼吸中止症評估方法及其系統/ METHOD ADAPTED TO DIAGNOSE AIRWAY OBSTRUCTION AND SYSTEM THEREOF	EU (Germany) Japan Taiwan USA
預測睡眠呼中止症風險之方法/ METHODS FOR PREDICTING THE RISK OF OBSTRUCTIVE SLEEP APNEA	USA Taiwan
[USS] 利用散射子分布統計量分析超音波回音信號之方法 ANALYSIS METHODS OF ULTRASOUND ECHO SIGNALS BASED ON STATISTICS OF SCATTERER DISTRIBUTIONS	Taiwan China USA EU (Germany)
[USV] 聲波散射訊號之影像強化方法與成像加速系統 [TW] IMAGE ENHANCEMENT METHOD AND IMAGING ACCELERATION SYSTEM FOR ULTRASOUND SCATTERER SIGNAL [US]Acceleration and enhancement methods and system for ultrasound scatterer structure visualization	USA
[UV] 超音波都卜勒影像之分群、雜訊抑制及視覺化方法 CLUSTERING, NOISE REDUCTION AND VISUALIZATION METHOD FOR ULTRASOUND DOPPLER IMAGES	Taiwan USA
細胞學影像特徵量化與視覺化方法及系統/ Method and system for quantification and visualization characteristics of cytology images	Taiwan USA

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
JAMA Otolaryngology - Head & Neck Surgery
 Association of Backscattered Ultrasonographic Imaging of the Tongue With Severity of Obstructive Sleep Apnea in Adults. vol. 149,7 (2023): 580-586.

Pattern Recognition
 Multivariate multi-layer classifier. vol. 131, 108896, July 2022.


Sleep Medicine
 Using standardized ultrasound imaging to correlate OSA severity with tongue morphology. vol. 120 15-21. 31 May. 2024.


Ultrasonics
 Clinical performance of ultrasonic backscatter parametric and nonparametric statistics in detecting early hepatic steatosis. vol. 142, 107391. Aug. 2024.


European Thyroid Journal
 Reliability of a computer-aided system in the evaluation of indeterminate ultrasound images of thyroid nodules. vol. 11,1. e210023. 1 Jan. 2022.


Ultrasound in Medicine & Biology
 Diagnostic performance of ultrasound computer-aided diagnosis software compared with that of radiologists with different levels of expertise for thyroid malignancy: A multicenter prospective study. vol. 47,1 (2021): 114-124.

Cancers
 Multi-Reader Multi-Case Study for Performance Evaluation of High-Risk Thyroid Ultrasound with Computer-Aided Detection. vol. 12,2 373. 6 Feb. 2020.

Clinical Endocrinology
 Differences in the ultrasonographic appearance of thyroid nodules after radiofrequency ablation. vol. 95,3 (2021): 489-497.

Original Research—General Otolaryngology

Backscattered Ultrasonographic Imaging of the Tongue and Outcome in Hypoglossal Nerve Stimulation

Samuel Tschopp, MD^{1,2*}, Vlado Janjic, MD^{1*}, Yili Lee, MBA³, Argon Chen, PhD⁴, Pei-Yu Chao, PhD³, Marco Caversaccio, MD², Urs Borner, MD^{2**}, and Kurt Tschopp, MD^{1**}

Abstract

Objective. Hypoglossal nerve stimulation (HNS) is an increasingly used therapy. However, not all patients undergoing HNS implantation benefit from the treatment, making an improved patient selection a priority. This study investigates whether backscattered ultrasonographic imaging (BUI) can predict the response to HNS therapy.

Study Design. Cross-sectional study.

Setting. Secondary and tertiary hospital.

Methods. In this multicenter cross-sectional study, we recruited patients who had undergone HNS implantation during their scheduled follow-up consultation. HNS therapy parameters were collected. Standardized submental ultrasonographic examination and home sleep apnea testing were performed. The primary outcome was assessing the response to HNS therapy using ultrasonographic features and preoperative patient characteristics.

Results. In total, 62 participants, 49 male, with a median (interquartile range [IQR]) age of 62 (55-67) and a median (IQR) body mass index of 27.6 (25.2-29.7). The follow-up was a median (IQR) of 19.5 (4.8-41.4) months after implantation. The apnea-hypopnea index (AHI) was preoperatively 40.5 (29.8-58.0) and reduced at follow-up to 21.0 (11.0-35.3). In total, 42% were responders to HNS. Preoperative AHI (34.8/hour vs 49.3/hour; $r = 0.44$) was significantly higher in nonresponders than in responders. The average prediction accuracy of HNS therapy based on baseline AHI alone was 71%. A lower backscatter signal, indicating less fat deposition in the tissue, was observed in the responder group. When the baseline AHI and backscatter signal were combined, the prediction accuracy of response to the HNS reached 78%.

Conclusion. The combination of tissue composition analyzed using the backscattered signal and the preoperative AHI is highly predictive for determining the HNS treatment response.

Trial Registration. ClinicalTrials.gov identifier NCT06154577.

Otolaryngology—
Head and Neck Surgery
2025, Vol. 00(00) 1–7
© 2025 The Author(s).
Otolaryngology—Head and Neck
Surgery published by Wiley
Periodicals LLC on behalf of
American Academy of
Otolaryngology—Head and Neck
Surgery Foundation.
DOI: 10.1002/oto.1251
http://otojournal.org
WILEY

Keywords

hypoglossal nerve stimulation, obstructive sleep apnea, outcome, ultrasonography

Received November 18, 2024; accepted March 23, 2025.

Hypoglossal nerve stimulation (HNS) is increasingly important in managing patients with obstructive sleep apnea (OSA) who do not tolerate continuous positive airway pressure therapy and are not eligible for other alternative treatment options, such as mandibular advancement devices or positional therapy.¹ The number of HNS implantations is increasing, with several thousand new implantations per year, and they have found their place in treatment guidelines.² Even with restrictive patient selection according to guidelines and optimal HNS titration, the responder rate remains around 50% to 60%.^{1,3-6}

Current anatomical and physiological selection criteria for respiration-triggered HNS include moderate to severe OSA, intolerance to continuous positive airway pressure therapy, body mass index (BMI), drug-induced sleep endoscopy (DISE), and exclusion of central sleep apnea,

¹Department of Otorhinolaryngology, Head and Neck Surgery, Kantonsspital Baselland Liestal, Liestal, Switzerland

²Department of Otorhinolaryngology, Head and Neck Surgery, Inselspital, University Hospital and University of Bern, Bern, Switzerland

³AmCad Biomed Corporation, Taipei, Taiwan

⁴Graduate Institute of Industrial Engineering, National Taiwan University, Taipei, Taiwan

*These authors contributed equally to this article and are equally contributing first authors.

**These authors contributed equally to this article and are equally contributing last authors.

Corresponding Author:

Samuel Tschopp, MD, Department of Otorhinolaryngology, Head and Neck Surgery, Inselspital, University Hospital and University of Bern, Freiburgrasse 20, 3010 Bern, Switzerland.
Email: samuel.tschopp@insel.ch

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