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

Scatterer Structure Visualization & Ultrasound Echo Analysis

Tumor Contour Retrieval

Multi-Layer Classifier & Echogenicity /Tumor Feature /Echo Texture Quantification

Other Patent Titles	Countries
用於頭頸部評估或介入治療的定位裝置/ POSITIONING APPARATUS FOR HEAD AND NECK ASSESSMENT OR INTERVENTION	China
	EU (Germany)
	Japan
	Taiwan
	USA
呼吸中止症評估方法及其系統/ METHOD ADAPTED TO DIAGNOSE AIRWAY OBSTRUCTION AND SYSTEM THEREOF	China
	EU (Germany)
	Japan
	Taiwan
	USA
COMPUTER-IMPLEMENTED METHOD FOR PREDICTING THE RISK OF OBSTRUCTIVE SLEEP APNEA	USA
	Taiwan
	EU
[USS] 利用散射子分布統計量分析超音波回音信號之方法 ANALYSIS METHODS OF ULTRASOUND ECHO SIGNALS BASED ON STATISTICS OF SCATTERER DISTRIBUTIONS	Taiwan
	China
	USA
	EU
	(Germany)
[USV] 聲波散射訊號之影像強化方法與成像加速系統 [TW] IMAGE ENHANCEMENT METHOD AND IMAGING ACCELERATION SYSTEM FOR ULTRASOUND SCATTERER SIGNAL	USA
[US]Acceleration and enhancement methods and system for ultrasound scatterer structure visualization	
[UV] 超音波都卜勒影像之分群、雜訊抑制及視覺化方法 CLUSTERING, NOISE REDUCTION AND VISUALIZATION METHOD FOR ULTRASOUND DOPPLER IMAGES	Taiwan
	USA
細胞學影像特徵量化與視覺化方法及系統/ Method and system for quantification and visualization characteristics of cytology images	Taiwan
	USA

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Association of Backscattered Ultrasonographic Imaging of the Tongue With Severity of Obstructive Sleep Apnea in Adults. "JAMA Otolaryngology Head & Neck Surgery" 2023, 149(7), 580-586.


Using Standardized Ultrasound Imaging to Correlate OSA Severity with Tongue Morphology. "Sleep Medicine" 2024, Volume 120, 15-21.


Diagnostic performance of ultrasound computer-aided diagnosis software ... "Ultrasound in Medicine and Biology" 47(1), 114-124.

Re-liability of a computer-aided system in the evaluation of indeterminate ultrasound images of thyroid nodules. "ETJ" 21-0023.


Multi-Reader Multi-Case Study for Performance Evaluation of High-Risk Thyroid Ultrasound....in "Cancers" 2022, 12(2), 373.


Ultrasound Backscatter Envelope Statistics Parametric Imaging for Liver Fibrosis in "Ultrasonic Imaging" Volum 42, Issue 2


Computer-aided diagnosis system of thyroid nodules ultrasonographyin "Medicine" 99(23), e20634.

Differences in the ultrasonographic appearance of thyroid nodules after radiofrequency ablation in "Clinical Endocrinology" Volume 95, Issue 3

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CDMO / Licensing

Use our expertise in developing product according to our strategic partners' needs/specifications.

Distribution

Leasing

Revenue
Sharing

Outright
Sales

Subscription

Exclusive or regional rights, single & multi-product deals.

Collaborative Marketing

Collaborate with enterprise imaging solutions providers to expand market access.

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Using standardized ultrasound imaging to correlate OSA severity with tongue morphology

Pien F.N. Bosschieter^a  , Stanley Y.C. Liu^a, Pei-Yu Chao^b , Argon Chen^c ,
Clete A. Kushida^d  

^a Department of Otolaryngology, head and Neck Surgery, Stanford University School of Medicine, Stanford, CA, USA

^b AmCad BioMed Corporation, Taipei, Taiwan

^c Institute of Industrial Engineering, National Taiwan University, Taipei, Taiwan

^d Department of Psychiatry & Behavioral Sciences, Stanford University School of Medicine, Stanford, CA, USA

Received 18 March 2024, Revised 21 May 2024, Accepted 28 May 2024, Available online 31 May 2024, Version of Record 5 June 2024.

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Collaborative Research and Supply Agreement

This Agreement is made and entered into as of Feb. 29th 2024 (the "Effective Date") by and between AmCad BioMed Corporation, a corporation with office at 5F., No.167, Fuxing N. Rd., Songshan Dist., Taipei City 105, Taiwan R.O.C (hereinafter referred to as "AmCad") and the University of the Pacific, Arthur A. Dugoni School of Dentistry located at 155 5th St, San Francisco, CA 94103, United States ("UOP"). AmCad and UOP are individually referred to as the "Party"; or collectively referred to as the "Parties".

1. **Product Definition:** The "Product" referred to in this Agreement encompasses the AmCAD-UO solution, which integrates both AmCAD-UO software and the specified accompanying hardware as detailed in the table below:

Item Description	Qty
AmCad-UO computer aided detection (CAD) software	1
Head & Neck Laser Alignment Robotic Arm system	1
Terason uSmart 3200t ultrasound model	1
5C2A Curved probe (for AmCAD-UO application)	1

2. **Scope of Product Use:** UOP agrees to use the Product for clinical, teaching, and academic research purposes. This includes incorporating the Product into clinical trials, educational programs, and research studies conducted by UOP.

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Meeting AmCAD-UT Software - Geneva Hospital

Fri, Jun 14, 2024 | 4:00pm - 4:21pm | [Premium report](#)

Tags (1)

Participants (7)

Download

Share

Recap Deep Dive Coaching

Feedback 0 0 5

Notes Transcript

Search Settings Copy Edit

Read Score 83 GOOD Engagement 87 GOOD Sentiment 79 OKAY Read Average

Summary


The participants discussed technical issues related to Teams and note-taking AI, including the integration of a software into their system and addressing security concerns. They also explored the format of the thyroid report, with a focus on automatic report generation and the process of sending patient information along with the report. The purchasing department requested an official offer for a software and expressed interest in testing it for two years, with the possibility of extending the quotation... [View more](#)

Chapters & Topics

Chapter descriptions

- Technical Integration and Security Concerns 0:21
 - Practical aspects of using the software in a hospital
 - Pricing and contract terms
 - Installation and security requirements

Discussion on Thyroid Report Format 0:00

Highlights 14 Chapters 3

Reaction 1:59 Jump to

Decrease in Engagement
"It's something from your company. Where is that? Oh, yeah, yeah, yeah. Well, that is ..."

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13.-14. Juni 2024 | 13-14 juin 2024

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Predictive Value of Standardized Submental Ultrasonography for Response to Hypoglossal Nerve Stimulation

Samuel Tschopp^{1,2}, Vlado Janjic², Yili Lee², Argon Chen³, Pei-Yu Chao³, Marco Caversaccio¹, Urs Borner², Kurt Tschopp²

1 Department of Otorhinolaryngology, Head and Neck Surgery, Inselspital, Bern
2 Department of Otorhinolaryngology, Head and Neck Surgery, Kantonsspital Baselland, Liestal
3 AmCAD Biomed, Taipei, Taiwan

Objective

The overall reported hypoglossal nerve stimulation (HNS) response rate for apnea-hypopnea index (AHI) reduction is approximately 60%. The present study aims to evaluate the predictive value of standardized submental ultrasonography.

Methods


We performed ultrasonography in patients with HNS and healthy controls using the AmCAD-UO device. Anatomy of the soft tissue anterior to the upper airway and the tissue morphology using backscatter ultrasound images were measured. Backscatter imaging relies on scattering of ultrasound waves by tissue inhomogeneities to quantify tissue composition and structure. Responder status for AHI-reduction was defined according to the Sher criteria.

Results


We enrolled 49 patients, of which 22 were responders and 27 were non-responders. In the HNS group, responders and non-responders were comparable in sex and age. However, the initial neck circumference (42.2 vs. 40.4cm), baseline AHI (50.8 vs. 34.7/h), and BMI (28.4 vs. 26.4kg/m²) were significantly higher in non-responders versus responders. The prediction of responders based on baseline AHI alone was 64%. A lower backscatter ultrasound image value, indicating less fat deposition in the tissue, was observed in the responder group. The backscatter ultrasound image from the oropharynx region, combined with baseline AHI, predicted responders with an accuracy of 81% and non-responders with 75% accuracy. When the baseline AHI and backscatter signal were combined, the area under the curve to predict treatment response to the HNS reached 0.80.

AmCAD ultrasonography with the robotic probe positioning for a standardized examination of the upper airway


B-mode ultrasonography of the oropharynx with various anatomical structures marked (left). For reference, an MRI image of the assessed airway levels during standardized ultrasonography

Qualitative assessment Quantitative assessment

Axial scans of a patient with obstructive sleep apnea and hypoglossal nerve stimulation in B-mode for qualitative assessment and backscattered signal for quantitative analysis.

Backscattered signal of the oropharynx region and baseline apnea-hypopnea index (AHI) predict responder to hypoglossal nerve stimulation

Conclusion

The combination of tissue morphology analyzed using standardized submental ultrasonography and patient characteristics is highly predictive for determining the treatment response in patients with HNS.

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Kantonsspital Baselland

Connect with Samuel Tschopp for more details or discussion:
samueltschopp@insel.ch

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
7 of 19

Article

Relationship between Cephalometric and Ultrasonic Airway Parameters in Adults with High Risk of Obstructive Sleep Apnea

Anutta Terawatpothong¹, Chidchanok Sessrisombat¹, Wish Banhiran², Hitoshi Hotokezaka³, Noriaki Yoshida³ and Irim Sirisoontorn^{1,4,*}

- ¹ Department of Clinical Dentistry, Walailak University International College of Dentistry (WUICD), 67 Ranong 2 Road, Dusit, Bangkok 10300, Thailand
- ² Department of Otorhinolaryngology, Faculty of Medicine Siriraj Hospital, Mahidol University 2 Wanglang Road, Bangkoknoi, Bangkok 10700, Thailand
- ³ Department of Orthodontics and Dental Facial Orthopedics, Nagasaki University Graduate School of Biomedical Sciences, 1-7-1 Sakamoto, Nagasaki 852-8588, Japan
- ⁴ Correspondence: irim.sirisoontorn@gmail.com

Abstract: Background/Objectives: Polysomnography and cephalometry have been used for studying obstructive sleep apnea (OSA) etiology. The association between craniofacial skeleton and OSA severity remains controversial. To study OSA's etiology, cephalometry, fiberoptic pharyngoscopy, polysomnography, and sleep endoscopy have been used; however, airway obstructions cannot be located. Recent research suggested ultrasonography for OSA screening and upper airway obstruction localization. Thus, this study aims to investigate the relationship between specific craniofacial cephalometric and ultrasonic airway parameters in adults at high risk of OSA. Methods: To assess craniofacial structure, lateral cephalograms were taken from thirty-three adults over 18 with a STOP-Bang questionnaire score of three or higher and a waist-to-height ratio (WHtR) of 0.5 or higher. Airway parameters were assessed through submental ultrasound. Results: NSBA correlated with tongue base airway width, while MP-H correlated with oropharynx, tongue base, and epiglottis airway width. SNA, SNB, and NSBA correlated with tongue width at the oropharynx. At tongue base, ANB and MP-H correlated with tongue width. SNB and NSBA were associated with deep tissue thickness at the oropharynx, while MP-H correlated with superficial tissue thickness at velum and oropharynx. Conclusions: Cephalometric parameters (SNA, SNB, ANB, NSBA, and MP-H) were correlated with ultrasonic parameters in the velum, oropharynx, tongue base, and epiglottis.

Keywords: obstructive sleep apnea; lateral cephalometry; ultrasonography

1. Introduction

Obstructive sleep apnea (OSA) is a breathing disorder caused by an upper airway constriction that prevents normal ventilation during sleep [1]. Respiratory effort increases to sustain airflow via a constricted airway, and this is followed by an increase in serum carbon dioxide (hypercapnia) and a drop in serum oxygen (hypoxemia). Cortical arousal is attributed to the increased work of breathing during sleep. Common symptoms of OSA include snoring, gasping, or choking respiration; frequent nighttime awakenings; excessive daytime sleepiness; and irritability. Many previous studies have reported that there are many factors that increase the risk of OSA, which include male sex, obesity (higher BMI), older age, greater neck and waist circumference, as well as craniofacial disposition and pharyngeal soft tissues [2–4]. Untreated obstructive sleep apnea (OSA) can lead to severe problems such as coronary artery disease, congestive heart failure, myocardial infarction, hypertension, stroke, cardiac arrhythmia, and mortality [5].

For diagnosing obstructive sleep apnea, polysomnography (PSG) is the most reliable technique [1,7]. The American Academy of Sleep Medicine (AASM) has proposed a classification system for the severity of OSA based on the apnea-hypopnea index (AHI).

OSA is classified as mild to moderate (AHI ≥ 5 events/h and AHI < 30 events/h) or severe (AHI ≥ 30 events/h). Without conducting home sleep apnea testing or PSG, the use of clinical devices, questionnaires, and prediction algorithms is not recommended to diagnose OSA in adults. Uncomplicated adult patients who present with signs and symptoms indicating an increased risk of moderate-to-severe OSA can be diagnosed using home sleep apnea testing with a suitable device or PSG, while adults with comorbid conditions are recommended to use PSG for OSA diagnosis rather than home sleep apnea testing [1]. Following diagnosis, a treatment plan should be developed, which may involve behavioral therapy, oral appliances, positive airway pressure devices, surgery, and/or additional treatments. OSA is regarded as a chronic condition that needs multidisciplinary therapy [1,6].

However, polysomnography (PSG) is costly, time-consuming, many sleep laboratories have long waiting lists of patients, and an unfamiliar environment may lead to discomfort and perhaps compromise sleep quality [7]. A previous study reported that PSG, fiberoptic pharyngoscopy, and sleep endoscopy have all been used to investigate the predisposing etiological variables of OSA. Nevertheless, even with the use of all available methods, there are limitations in accurately locating the site of the blocked airway because it is not possible to mimic the sleep state in awake individuals [8]. PSG gives a range of physiological data during sleep without determining the anatomic cause of OSA. Consequently, a significant number of people remain undiagnosed due to an inability to access treatment, being asymptomatic, or lacking awareness of their condition. It is essential to place importance on the screening procedures that help in the early identification of patients at risk or those with unrecognized OSA. The screening methods for OSA involve physical examination, upper airway examination, questionnaires, and imaging technologies.

The STOP-Bang questionnaire, Berlin Questionnaire (BQ), and Epworth Sleepiness Scale (ESS) are widely used screening tools facilitating the early identification of patients at risk [9]. The STOP-Bang questionnaire is a reliable tool for screening obstructive sleep apnea in adults. It consists of four subjective items (STOP: Snoring, Tiredness, Observed Apnea, and High Blood Pressure) and four demographic criteria (Bang: BMI > 35 kg/m², age > 50 years, neck circumference > 40 cm, and gender). With a high sensitivity, the STOP-Bang questionnaire can identify patients with moderate-to-severe OSA [10,11]. However, some adjustments are required for Thai patients. The Thai STOP-Bang questionnaire's optimal BMI cutoff point was > 30 kg/m², which increased sensitivity at the AHI cutoff points of 5 and 15 by 88.7% and 93.2%, respectively. A patient is considered high-risk if at least three "yes" answers are provided [12]. Banhiran et al. reported that male gender and a waist-to-height ratio (WHtR) of 0.55 were significant predictors of moderate-to-severe OSA. The WHtR identifies central obesity, or visceral fat, and represents fat distribution [13]. However, some authors reported that if the WHtR is higher than 0.5 or 0.6, the risk of causing obstructive sleep apnea increases [14–17]. Therefore, the waist-to-height ratio (WHtR) should be incorporated into questionnaires for high-risk individuals to reduce the possibility of false-positive results.

Imaging technologies, such as cone-beam computed tomography (CBCT), lateral cephalometric radiographs, and ultrasonography, aid in the study of OSA. Cephalometric radiographs are routinely taken for orthodontic diagnosis and treatment planning by analyzing the relationship between the skull, maxilla, mandible, teeth, and soft tissues. The radiation dose from a lateral cephalogram is low, at approximately 5.6 μ Sv, whereas the CBCT is 15–26 times higher [18]. Previous research has found that most common variables, such as the angle between the anterior cranial base and the maxilla (SNA), the angle between the anterior cranial base and the mandible (SNB), the posterior airway space (PAS), the length of soft palate (PNS-P), and the distance from the hyoid bone to the mandibular plane (MP-H), may be correlated with the development and severity of OSA [19]. In addition, Banhiran et al. reported that the risk of moderate-to-severe OSA in Thai patients was also influenced by the cranial base angle (NSBA), which was 130 degrees or less. Abnormality of the skull's base may potentially contribute to the constriction of


Figure 4. Submental ultrasonography equipment, laser alignment (AmCad BioMed Corporation).

The ultrasonic transducer is aligned with the HM plane's laser projection in the submental region to provide transverse, cross-sectional images (Figure 4). Three replicates were automatically scanned by an ultrasonic transducer during tidal breathing and Muller maneuvers at 30 degrees, covering the velum (V), the oropharynx (O), the base of the tongue (T), and the epiglottis (E) (Figure 5a).


Figure 5. (a) The 30° segment of the upper airway; HM, hyoid-external meatus. (b) Transverse view of ultrasonographic images.

During scanning, patients were instructed to alternate between regular breathing and Muller maneuvers, completing a total of three repetitions. Muller maneuvers were conducted by inhaling vigorously with the mouth and nose close, simulating upper airway obstruction during awake.

The airway's width, airway's depth, tongue width, airway-tissue width ratio, superficial tissue thickness, and deep tissue thickness were measured (Figure 5b). Dynamic changes in airway width during tidal breathing and Muller maneuvers, or percentages of airway contraction, were calculated as a percentage of pharyngeal diameters shortening from the diameters during tidal breathing to that during MM.

The risk of moderate-to-severe OSA (AHI ≥ 15 /h) was evaluated as a percentage based on the assessment of ultrasonic parameters. The risk of 0–40% is classified as low, 40–80% is classified as medium, and 80–100% is classified as a high risk of moderate-to-severe OSA [24,27–29].

Citation: Terawatpothong, A.; Sessrisombat, C.; Banhiran, W.; Hotokezaka, H.; Yoshida, N.; Sirisoontorn, I. Relationship between Cephalometric and Ultrasonic Airway Parameters in Adults with High Risk of Obstructive Sleep Apnea. *J. Clin. Med.* 2024, 13, 3540. <https://doi.org/10.3390/jcm13123540>

Academic Editor: Jari P. Aikberg


Received: 16 May 2024
Revised: 11 June 2024
Accepted: 11 June 2024
Published: 17 June 2024

Copyright: © 2024 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<https://creativecommons.org/licenses/by/4.0/>).

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Department of Otolaryngology

Walailak University

International College of Dentistry
Clinical Dentistry Program

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
AI for POCUS

AI-powered tools on Venue™ family ultrasound systems help increase clinician confidence by detecting and tracking nerves during scouting for nerve block procedures, while Caption Guidance helps facilitate acquiring cardiac images in real time.

[Learn more >](#)

AI for Primary Care

Whizz Label™ on Versana Premier™ and Versana Balance™ automatically labels liver, gallbladder, and right kidney on images during abdomen scans of the RUQ (right upper quadrant), helping to save time and enhance workflow productivity.

[Learn more >](#)

AI for General Imaging

LOGIQ™ ultrasound systems now feature Auto Abdominal Suite 1.0. These three new onboard AI-powered, deep learning-based algorithms are designed to reduce intra-/inter-operator variability, enhance exam consistency, and increase workflow efficiency.

[Learn more >](#)

AI for Cardiology

Vivid™ devices offer AI-supported tools like Easy AFV LV for near-instant strain measurement and Easy AutoEF for one-click ejection fraction assessment, in AI-aided review for urgent heart issues for precision, speed and ease of use measurement.


[Learn more >](#)

AI for Women's Health

Reduce 2nd trimester exam time by up to 40% with SonoLyst™, a suite of tools available on Voluson™ systems, that leverage AI to identify fetal anatomy seen on standard views, then automatically annotates and measures where applicable.

[Learn more >](#)

AI for Breast Ultrasound

Breast Assistant, powered by Kolos D5™ on the LOGIQ E10 and Invenia™ ABUS 2.0 ultrasound systems, automatically populates a quantitative breast malignancy risk assessment aligned to BI-RADS® ATLAS category in as little as two seconds.

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摘自Key Market Trends in Handheld Ultrasound

Handheld **ultrasound** has not yet reached mainstream adoption, but the market is still forecast to reach over \$500 million by 2026. Signify Research's newly-published **Handheld Ultrasound Deep Dive Report 2022** shows the market revenues are estimated to have grown by approximately 30% in 2021, driven by strong growth in the U.S., the biggest market for handheld ultrasound. Despite the global challenges for handheld ultrasound vendors in 2022, such as rising inflation and supply chain disruptions, the handheld ultrasound market is expected to experience double-digit growth, and this is forecast to continue through to 2026. Most of the market growth will be fueled by increased adoption of handheld devices by new users of ultrasound, such as primary care physicians, nurses, emergency medical technicians (EMTs) and midwives. The key market trends are discussed below.

Trend 1: Artificial Intelligence

The steep learning curve and subsequent ultrasound skills shortage are two of the biggest barriers to the wider use of ultrasound. These challenges are exacerbated in **handheld ultrasound**, where there is a higher proportion of new and less experienced users compared to cart and compact ultrasound. As the expansion of the handheld ultrasound market in the coming years is expected to be strongly driven by new user groups, this barrier will be greater than ever before. These barriers can be partially addressed by **artificial intelligence** (AI) solutions that guide users with positioning and moving the ultrasound probe.

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